

COVID19 Rapid Test (nasal) for pre-school children

Step 1: Preparation

- Wash or disinfect hands
- Get material ready



Step 2: Prepare material

- Open test kit
- Get tube with buffer solution and nozzle cap ready
- Remove swab from the package
- Get paper tissue ready


Step 3: Nasal swab

- Insert swab ca. 1-2cm
- Rotate swab 4 times (ca. 5 sec.) at the nasal wall
- Repeat in the other nostril




Step 4: Swab into tube

- Take the tube into the other hand and take off the cap
- Put the swab into the solution
- Squeeze the tube
- Rotate the swab ca. 10 times
- Remove the swab while squeezing it and place it onto the paper tissue




Step 5: Drop onto test strip

- Place the nozzle cap onto the tube
- Apply 4 drops onto the test strip, wait for 15 min.



Step 6: Interpret result

C= Control line
T= Test line



Test is **negative** = no indication of an infection with SARS-CoV2

Test is **positive** = Infection with SARS-CoV2

Test is **invalid** = Repetition immediately required

Please dispose of all materials in the prepared waste bins and then disinfect or wash hands!